

Client Notification

Date: August 30, 2021

From: Wade K. Aldous, PhD, Disease Control Director, State Hygienic Laboratory

Effective Date: October 1, 2021

Subject: SARS-CoV-2 IgG Antibody testing changes

Reason: Manufacturer has developed a new semi-quantitative reagent as their next generation in COVID-19 serology testing. The “old” reagent will no longer be manufactured.

DiaSorin Old test result interpretation:

- Negative < 15.0 AU/mL
- Positive ≥ 15.0 AU/mL

DiaSorin New test result interpretation:

- Negative < 13.0 AU/mL
- Positive ≥ 13.0 AU/mL

Thank you,

Wade K. Aldous, PhD, Disease Control Director

Wade-aldous@uiowa.edu

319-335-4765

Molly Bradshaw, Clinical Lab Technical Specialist

Molly-bradshaw@uiowa.edu

319-335-4291